Analysis of Soils and Sediments for Loss on Ignition

USGS-Mercury Research Laboratory 8505 Research Way Middleton, Wisconsin 53562 <u>mercury@usgs.gov</u> (608) 821-3844

Method Summary

The following standard operating procedure (SOP) describes the analysis of soils and sediments for loss on ignition (LOI) by the U.S. Geological Survey's Mercury Research Laboratory (MRL). Sample is weighed into aluminum boats and heated to 550°C for two hours. The percent of sample mass lost following heating is reported as LOI. One sample out of every 14 analyzed is weighed in triplicate to assess method precision.

Laboratory Safety

Analysts who use the MRL must have read, understood, and signed the Chemical Hygiene Plan for the MRL prior to potential exposure to any chemicals. The analyst must have a thorough understanding of the required safety protocols for the lab chemicals prior to their use of the lab. Adequate personal protection equipment such as safety glasses, gloves, and chemical resistant clothing must be worn when exposure to hazardous chemicals are possible. Caution should always be exercised; chemicals are present in the laboratory and often in use by other analysts. Hazardous chemicals should only be handled by adequately trained personnel under a high volume fume hood with extreme caution.

The analysis of LOI involves extremely high temperatures, and due caution should be exercised. Adequate personal protection equipment such as safety glasses, heat shielding gloves, and chemical resistant clothing must be worn when required.

Sample Analysis

Unless otherwise specified, solid samples are lyophilized ("freeze-dried" under vacuum while frozen) to a consistent weight and homogenized (via ball mill, coffee grinder/food processor, or mortar/pestle) prior to analysis. The samples should appear to be well pulverized and mixed to a consistent composition before weighing. A typical analysis contains up to 64 samples, with one sample out of every 14 weighed in triplicate

- 1. Arrange an adequate number of aluminum boats onto a tray. The maximum number of boats the furnace can accommodate is 64.
- 2. Using a ball point pen, number each aluminum boat on the tag. Make sure that you press hard enough to leave an imprint; the ink will not be visible following combustion.
- 3. Weigh 100 500 mg of sample into each aluminum boat. Make sure to record into the appropriate data sheet the vial ID, aluminum tray number, empty aluminum boat mass, and the mass of the aluminum boat with the added sample.

- Carefully transfer the boats to the furnace, place in the bottom, and heat at 550°C for two hours. During the combustion, open the adjacent lab door several inches and turn on the fume hood to provide ventilation out of the building.
- 5. After two hours, turn off the furnace and allow the samples to cool. Carefully remove the samples.
- 6. Weigh the combusted sample and record the combusted weight into the data sheet.

Quality Assurance and Control Objectives

A successful run must meet the following criteria. The relative standard deviation of the triplicate analysis must be < 10%, and all values for LOI must be positive. If the analysis fails either of these criteria, check that the combusted sample mass was entered into the correct position in the data sheet and/or reweigh the samples (if possible). If this does not correct the failure, repeat the entire analytical batch as that the data is likely compromised.

APPENDIX 1. Example of a completed LOI sample setup sheet.

Sample ID		Heade informa	tion Ta	Tin Tare Weight I			Sample ight Tin + Baked Sample Weight			
			Creek Seds.xls rmat <u>T</u> ools <u>D</u> ata	Window	<u>H</u> elp			_		
: 0		100			 8. 2↓ X↓	100%	• 🕜 🔡 🗄	Arial	- 10	
	H70 ·	∱ <mark>9</mark> RS								
	A	B	C	D	E		F (
1	Sample Set:	Bonanza Creek	Seds				Date: 2/23	3/2015		
2	Analy Analy	st: CD1		↓	1		V			
					Tin		mple			
4	Vial ID	Sample	D Tin ID	Tin W	eight Wei	-	aked eight % I	0		
	MSC930X	Sample	1	0.42	-		-	23308 % RSD		
5 6	MSC930X		2	0.42				76829 0.19%		
7	MSC930X		3	0.41				80263	,	
8	MSC925X		4	0.42				27027		
9	MSC934X		5	0.41				75723		
	MSC922X		6	0.41				56071		
	MSC932X		7	0.42			4298 0.875			
	MSC914X		8	0.42				46193		
	MSC917X		9	0.41				06935		
	MSC744X		10	0.42				54386		
	MSC943X		11	0.41				52708		
	MSC910X		12	0.41			419 0.964			
	MSC950X		13	0.42				71383		
	MSC939X		14	0.42				12844		
	MSC906X		15	0.41				54275		
	MSC951X		16	0.4				92095		
	MSC745X		17	0.41				60963 % RSD		
	MSC745X		18	0.42				60819 0.11%		
	MSC745X		19	0.41				74419		
	MSC908X		20	0.41			4208 0.977			
	MSC947X		21	0.42				18593		
	MSC928X		22	0.41			423 0.887			
	MSC748X		23	0.42			4226 0.973			
	MSC927X		24	0.41				68595		
	MSC902X		25	0.41				88889		
	MSC746X		26	0.42				34247		
	MSC944X		27	0.42			4225	0.915		
	MSC929X		28	0.42				84906		
	MSC750X		29	0.42				58824		
	MSC933X		30	0.41				25527		
	MSC936X		31	0.4				54098		
26	MSC001V	11 (Cha 12 (5)	20	0.42				02750		
Read		et1 / Sheet2 / She	eet3 /						•	

3

Appendix 2. Definition of equations.

$$\text{Loss on Ignition} = \frac{\left(\begin{pmatrix} \text{Mass of Tin} \\ \text{with Sample} \end{pmatrix} - \begin{pmatrix} \text{Tin Tare} \\ \text{Weight} \end{pmatrix} \right) - \left(\begin{pmatrix} \text{Mass of Tin} \\ \text{with Baked Sample} \end{pmatrix} - \begin{pmatrix} \text{Tin Tare} \\ \text{Weight} \end{pmatrix} \right)}{\left(\begin{array}{c} \text{Mass of Tin} \\ \text{with Sample} \end{pmatrix} - \begin{pmatrix} \text{Tin Tare} \\ \text{Weight} \end{pmatrix} \right)}$$